

Lipid Nanoparticle Production for mRNA Delivery: Comparison of Different Turbulent Mixing Technologies

Vogler J.

leon-nanodrugs GmbH, Munich, Germany

WEBSITE: www.leon-nanodrugs.com
EMAIL: j.vogler@leon-nanodrugs.com



PURPOSE

mRNA-loaded lipid nanoparticles (LNPs) have revolutionized vaccines & gene therapy drug products. Turbulent mixing has proven to be the preferred manufacturing technology for high-volume production due to its ease of scale-up compared to e.g. microfluidics.

Here, we show **comparative data on jet impinging technology from LEON vs. conventional T-junction mixing** to evaluate the impact of the mixing technology on the particle properties and *in vivo* transfection efficiency. A clinically relevant lipid composition loaded with polyA surrogate or firefly luciferase (FLuc) mRNA as reporter system was used.

METHODS

Mixing technologies:




Figure 1. LEON impinging jet mixer




Figure 2. Conventional T-junction mixer

- LNP production:** LNPs were produced by mixing a lipid solution consisting of ionizable cationic lipid/Cholesterol/DSPC/PEG-Lipid in ethanol with a buffered aqueous solution containing polyA or FLuc mRNA. A flow rate ratio (FRR) of 2:1, 3:1 and 4:1 as well as a total flow rate (TFR) of 40 ml/min and 80 ml/min was selected for the manufacturing of polyA loaded LNPs using the LEON impinging jet mixer. Results were compared with LNPs produced with a FRR of 3:1 and TFR of 40 ml/min using a conventional T-junction mixer. Subsequently Fluc mRNA LNPs were manufactured at a FRR of 3:1 and TFR of 40 ml/min with both LEON impinging jet mixer and conventional T-junction mixer. All samples were dialyzed against PBS pH 7.4 without prior dilution.
- Particle characterization:** Particle size & poly dispersity index (PDI) were determined by dynamic light scattering (DLS) using a Malvern Zetasizer (Zetasizer Nano ZS, Malvern, UK).
- Encapsulation efficiency:** Encapsulated mRNA was quantified using a fluorogenic RiboGreen™ RNA assay kit.
- In vivo transfection assay:** 5 mice each received a single IV (tail vein) injection of 0.3 mg/kg or 1 mg/kg mRNA-loaded LNPs produced with the LEON impinging jet reactor or a conventional T-junction mixer. Bioluminescence of liver homogenates was determined 4 h after LNP administration.

RESULTS

- Both mixing technologies generated polyA and FLuc-mRNA loaded LNPs with **comparable size, PDI, encapsulation efficiency (EE%)** (Table 1), **morphology** (Figure 3) and **in vivo activity** (Figure 2).




Figure 1. Size (z-average) and PDI of polyA and FLuc mRNA loaded LNPs measured by DLS post dialysis. Particles were produced at FRR=3:1 and TFR=40 ml/min with both LEON impinging jet reactor and conventional T-piece mixer.

Table 1. Size (z-average), PDI and encapsulation efficiency (EE%) of polyA and FLuc mRNA loaded LNPs produced with different process parameters. Results for particles produced with a conventional T-piece are highlighted in pink.

Sample	Process Step	Size (Z-Ave nm)	PDI	EE%	TFR (ml/min)	FRR	DS
Sample A (Impinging Jet)	Post-dialysis	71.6	0.093	95	40	2:1	Poly-A
Sample B (Impinging Jet)	Post-dialysis	61.78	0.082	97	40	3:1	
Sample C (Impinging Jet)	Post-dialysis	71.97	0.157	97	40	4:1	
Sample D (Impinging Jet)	Post-dialysis	62.01	0.103	98	80	3:1	
Sample E (T-piece)	Post-dialysis	72.33	0.091	99	40	3:1	FLuc
Sample F (Impinging Jet)	Post-dialysis	63.35	0.061	98	40	3:1	
Sample G (T-piece)	Post-dialysis	63.65	0.033	ND	40	3:1	
Sample H (Impinging Jet)	Post-conc.	65.68	0.007	99	40	3:1	
Sample I (T-piece)	Post-conc.	64.23	0.034	99	40	3:1	




Figure 2. In vivo transfection assay. Bioluminescence of liver homogenates was determined 4 h after IV injection with FLuc mRNA loaded LNPs produced with LEON impinging jet reactor or a conventional T-piece. n=5

- The variation of FRR for polyA loaded LNPs produced with LEON impinging jet reactor showed little impact on the particle properties with FRR=3:1 resulting in the smallest size and PDI. Increasing the TFR from 40 ml/min to 80 ml/min did not result in a substantial change in size and PDI when using the LEON impinging jet reactor.
- Both mixing technologies result in homogenous spherical particles according to cryo-TEM analysis.

CONCLUSIONS

- ✓ LEON's impinging jet reactor technology enables the manufacturing of mRNA loaded LNPs with **small size, low PDI and high encapsulation efficiency**.
- ✓ LNPs produced with LEON's impinging jet reactor technology are **comparable** with LNPs produced with a conventional T-piece mixer in terms of **particle size, PDI, EE%, morphology, and in vivo activity**.
- ✓ LNP characteristics (size, PDI, EE%) of particles produced with LEON's impinging jet reactor technology are **consistent over a wide range of total flow rates** (40 – 80 ml/min tested).
- ✓ The collected **scientific data demonstrate** once again the **quality** of LEON's process technology to support **product design and development** based on LEON's expertise and its **data driven approach**.

ACKNOWLEDGMENT

We thank Acuitas Therapeutics (Vancouver) for the collaborative support on the present work.

