

# Lipid Nanoparticle Production for mRNA Delivery: Impact of Mixing Technologies on Transfection Performance *in vivo*

Dubiella C.<sup>1</sup>, Soliman I.<sup>2</sup>, Tafa O. B.<sup>3</sup>, Balazs K.<sup>2</sup>, Gold L.<sup>2</sup>, Fois C.<sup>2</sup>, Rottensteiner H.P.<sup>2</sup>, Feichtinger G.A.<sup>2</sup>

<sup>1</sup>leon-nanodrugs GmbH, Munich, Germany

<sup>2</sup>Evotec GT GmbH, Orth an der Donau, Austria. <sup>3</sup>Aptuit, an Evotec Company, Verona, Italy.

WEBSITE: [www.leon-nanodrugs.com](http://www.leon-nanodrugs.com)  
EMAIL: [info@leon-nanodrugs.com](mailto:info@leon-nanodrugs.com)



## PURPOSE

mRNA-loaded lipid nanoparticles (LNPs) have revolutionized vaccines & gene therapy drug products. Two different mixing technologies exist for their production: **jet impinging** and **microfluidics**, whereas the former is preferred in **high-volume production** due to its ease of **scale-up**.

Here, we showcase a **side-by-side comparison** between **jet impinging (LEON equipment)** & **microfluidic technology** to evaluate the impact of the mixing principle on transfection performance **in vitro** & **in vivo** experiments. A clinically relevant lipid composition loaded with firefly luciferase (FLuc) mRNA as reporter system was used.

## METHODS

### • LEON Equipment train:



Equipment	Details	Volume output [ml/min]	Footprint (H x W x D) [m]
LEON 1	batchwise production	Bench scale: 1-80	1 x 1 x 0.5
LEON 2	continuous flow ready	Bench scale: 1-80	1 x 1 x 0.5
LEON 3	continuous flow ready	Pilot scale: 30-500	1 x 1 x 0.6
LEON GMP unit	continuous flow ready	Commercial scale: 300+	2 x 0.8 x 1.2

• **Formulation:** Onpatro® lipid composition: DLin-MC3/Cholesterol/DSPC/DMG-PEG2000 – 50/38.5/10/1.5 mol% with total lipid concentration of 10 mg/ml. FLuc mRNA (Trilink or APEX BIO) concentration of 0.09 mg/ml in non-solvent.

• **LNP production:** LNPs were produced with a flow rate ratio (FRR) of 3:1 (non-solvent:solvent) and total flow rate (TFR) = 30 ml/min using an impinging jet reactor setup (LEON 1) and TFR 10 ml/min using a microfluidics setup (Ignite™, Precision NanoSystems). Lipids were dissolved in EtOH (solvent) and precipitated against 50 mM citrate buffer (pH 4) (non-solvent). Samples were diluted, subjected to dialysis overnight (PBS pH 7.4) and concentrated to approx. 1 mg/ml.

• **Particle characterization:** Particle size & polydispersity index (PDI) were determined by dynamic light scattering (DLS) using a Stunner instrument (Unchained Labs). Cryo-transmission electron microscope (TEM) samples were inspected on a Tecnai F20 TEM.

• **Encapsulation efficiency:** Encapsulated mRNA was quantified using a fluorogenic Quant-iT™ RiboGreen™ RNA assay kit (Thermo Fisher Scientific).

• **In vitro transfection assay:** HepG2 cells were incubated in a 24-well plate format with various doses of FLuc mRNA-loaded LNPs.

• **In vivo transfection assay:** B6 albino mice (n=6 per group) received a single IV (tail vein) injection of 60 µg mRNA-loaded LNPs (2.4 mg/kg). Bioluminescence imaging on an IVIS Spectrum imaging system (PerkinElmer) was performed 6 h, 24 h & 48 h after LNP administration.

## RESULTS

- Both mixing technologies generated FLuc-mRNA loaded LNPs according to cryo-TEM (Figure 1) with **comparable size, PDI** (Figure 2) and **encapsulation efficiency (EE%)** (Figure 3).




Figure 1. Cryo-TEM image of a sample of FLuc mRNA-loaded LNPs produced with LEON's impinging jet technology showing spherical shaped particles.




Figure 2. Particle size (Z-average) & polydispersity index (PDI) of LNP samples analyzed by dynamic light scattering (DLS).




Figure 3. Encapsulation efficiency (EE%) of LNP samples determined by a fluorometric RiboGreen assay.




Figure 4. In vitro transfection assay using 10 µg mRNA (according to EE% assay) per transfection. All data n=6, non-parametric, Kruskal-Wallis post AVG+/-SD, \*\* p<0.01.




Figure 5. In vivo gene delivery as observed by luciferase bioluminescence imaging.




Figure 6. In vivo gene delivery quantification for total body and liver (top) and total body (bottom) over time.

- The high encapsulation efficiency (**EE% > 97%**) of samples of both mixing technologies translated into high transfection rates in cell assays (Figure 4) and *in vivo* (Figure 5 & 6). FLuc expression seemed to peak at 6h post dosing and declined thereafter. In vivo imaging suggests **systemic delivery of LNPs** in both groups (LEON 1 and microfluidics) and predominant transfection of the liver as target organ. **5/6 animals positive** in both groups.

## CONCLUSIONS

✓ **LEON equipment efficiently** enables the production of mRNA-loaded LNPs with **low PDI** (PDI < 0.12) and **high encapsulation efficiency** (EE% > 97%).

✓ LEON equipment enables the production of mRNA-loaded LNPs under low endotoxin burden RNase contamination.

✓ LNPs produced with LEON 1 equipment **under turbulent mixing conditions** (TFR: 30 ml/min) were **comparable** to LNPs produced with a **laminar flow-based system** (Ignite, Precision NanoSystems) in terms of:

- **Size & PDI** (72-79 nm, PDI < 0.12)

- **Encapsulation efficiency** (EE% > 97%).

- **In vitro & in vivo transfection established** - non-inferiority of LEON 1 vs. microfluidics produced mRNA/LNPs demonstrated based on expression yields and kinetics.

- Further studies with higher flow rates (TFR > 300 ml/min) are currently ongoing for high volume manufacturing of LNPs as part of LEON's database.

